



washroom cleaner probiotic

Ecological probiotic bathroom cleaner with long-lasting self sustaining effects. Especially suitable for hard-to-reach areas, such as (micro)pores, corners, and cracks. Regular use removes (gray) deposits, prevents dirt buildup, and combats unpleasant odors.

Better

Smart dosing options (automatic dispenser dosing, portable manual dosing, and ready-to-use) help prevent overdosing, reduce waste, and lower overall cleaning costs.

Safer

Non-classified and safe to use without protective equipment. Where industry or site-specific safety guidelines recommend the use of protective equipment, i-hygienic advises adhering to those recommendations.

Greener

Plant-based formula and bio-based packaging result in a net carbon-negative impact.

Faster

Performance-enhanced probiotic daily cleaner for quick and efficient washroom cleaning.

Cleaner

The optimal daily cleaning solution with probiotic deep-cleaning performance for washrooms and sanitary areas.

i.31 washroom cleaner probiotic

Instructions for Use

Prepare the cleaning solution according to the dosage instructions, using either a bucket or a spray bottle.

Bucket method: The prepared cleaning solution can be applied directly onto the surface. Allow it to work for 1–5 minutes if needed. For stubborn dirt, apply generously and loosen with a sponge or brush. Wipe off with a clean, damp cloth and let the surface air dry. In case of heavy soiling, the recommended dosage may be doubled. Spray method: Apply the prepared cleaning solution onto a damp cloth and use it to clean the surface. Wipe with a clean, damp cloth and allow to air dry. For floors, mop, wipe, or scrub with the cleaning solution. There is no need to rinse. Allow the floor to air dry.

Transport, Handling, and Storage

Non-classified product. No restrictions on handling, transport, or storage. See Safety Data Sheet (SDS) for detailed information.

Dosage

Dosing protocol	AutoDose Ultra	FlexDose Ultra	FlexDose	EasyDose
Packaging	Pouch	Dispensing bottle	Dispensing bottle	Spray bottle
Standard min. dosing	5ml on 5ltr	5ml on 5ltr	10ml on 5ltr	Ready to use*

*The product is ready to use; diluting it will reduce its cleaning effectiveness.

Specifications

Dosing protocol	AutoDose Ultra	FlexDose Ultra	FlexDose	EasyDose
Packaging (CU)	Pouch	Dispensing bottle	Dispensing bottle	Spray bottle
Content CU (ltr)	1.8	2	1	0.75
Articlecode CU	I31.PO.1800	I31.DB.2000	I31.DB.1000	I31.SB.750
Barcode CU	8721284400653	8721284400240	8721284400424	8721284401025
Articlecode box (BU)	K.4.I31.PO.1800	K.6.I31.DB.1000	K.6.I31.DB.1000	K.6.I31.SB.750
Barcode BU	8721284400752	8721284400325	8721284400523	8721284401285
Content BU	4	6	6	6
BU per layer	12	6	14	15
Pallet layers	5	5	6	5
BU on pallet	60	30	84	75
CU on pallet	240	180	504	450
Safety Data Sheet	Download SDS	Download SDS	Download SDS	Download SDS