


# iD.7 floor cleaner ultra degreasing



## Properties

Ecological ultra-concentrated degreasing floor cleaner for intensive machine cleaning of all water-resistant hard floor surfaces. Exclusively developed for use with the i-mop. Certified with the EU Ecolabel: DK/020/021.

### **Better**

Smart dosing options, pre-dosed pods or 2 L bottles with dosing cap, help prevent overdosing, reduce waste, and lower overall cleaning costs.

### **Safer**

Non-classified and safe to use without protective equipment. Where industry or site-specific safety guidelines recommend the use of protective equipment, i-hygienic advises adhering to those recommendations.

### **Greener**

Plant-based formula and biobased bottle result in a net carbon-negative impact.

### **Faster**

Performance-enhanced i-mop solution, offered as pre-dosed pods or in 2 L bottles with an integrated dosing cap for fast, consistent preparation.

### **Cleaner**

Strong degreasing performance removes oily residues and stubborn dirt from hard floors.

## product data sheet

# iD.7 floor cleaner ultra degreasing

### Instructions for Use

This product is intended for use with the i-mop only. Dose 10 ml of product using one dosing from the 2 L dosing bottle or one 10 ml i-dose pod directly into the i-mop tank. After dosing, shake gently to mix. Follow the manufacturer's guidelines for optimal and balanced results. After use of this cleaning product in the i-mop, drain the water tank, refill with fresh water, and clean the filters.

### Transport, Handling, and Storage

Non-classified product. No restrictions on handling, transport, or storage. See Safety Data Sheet (SDS) for detailed information.

### Dosage

Dosing protocol	FlexDose Ultra	i-dose
Packaging	Dispensing bottle	Pod
Standard min. dosing	5ml on 5ltr*	1 pod (10ml) on max. 5ltr*

\* The advised minimum dose applies to i-mop, i-fibre and similar high quality cleaning equipment; traditional methods may require a higher dose.

### Specifications

Dosing protocol	FlexDose Ultra	i-dose	i-dosemix (iD.3/iD.5/iD.7)
Packaging (CU)	Dispensing bottle	Pod	Pod
Content CU (ltr)	2	0,1	0,1
Articlecode CU	ID7.DB.2000	-	-
Barcode CU	8721284402596	-	-
Articlecode box (BU)	K.6.ID7.DB.2000	K.120.72.0444.6	K.120.72.0444.0
Barcode BU	8721284402633	8719778011799	8719778011973
Content BU	6	120	120
BU per layer	6	14	14
Pallet layers	5	10	10
BU on pallet	30	140	140
CU on pallet	180	16800	16800

**Safety Data Sheet**

[Download SDS](#)